


















Distributed by Cardinal Health, Sleep Products Division 800-582-7458
 800-742-3646 US and Canada Phone 603-758-6600 Fax 603-758-6699 www.sleep-net.com
 UK Office Phone/Fax +(44)186-532-7196

SleepNet Corp. MiniMe Headgear Sizing Chart

-  Recommended
-  May also fit



To find head size, measure circumference of head, just above eyebrows. Select the most appropriate size.

Head Size - Metric	35.5	38	41	43	46	48.5	49.5	52
Inches	14	15	16	17	18	19	19 1/2	20 1/4
								
								
								

MiniMe Large Mask and Headgear Part #TMS-8002
 Headgear Part #TMS-89MM

MiniMe Medium Mask and Headgear Part #TMS-8003
 Headgear Part #TMS-89MMPET

MiniMe Small Mask and Headgear Part #TMS-8004
 Headgear Part #TMS-89MMPXS